

User Manual



Company: Lifelab BV

✉ E-mail: sales@lifelab.shop

📍 Address: Kerstraat 107a, 5411CK Zeeland, Netherlands

for Facial

Instruction for Facial

Dear respected customer,

Thanks for choosing this product. We sincerely hope you are satisfied with our products and services. Please enjoy your lighting therapy according to the instruction below.

Please wear glasses or close eyes when enjoying therapy. Thank you very much.

Instruction of keys



Press 3s to turn on/off the device

Select mode
To change one mode to another, please pause first

Timing
(5/10/15/20mins)

Pause/Play

Rejuvenation

Rejuvenation

- Spectrum: Red 630nm&660nm
- Timing: 5/10/15/20 minutes
- Usage: 1-2 times a day, 5-30cm distance

Acne

Acne

- Spectrum: Blue 450nm/Red 630nm & 660nm/NIR 850nm&940nm
- Timing: 15minutes
- (5mins Blue + 10mins Blue/Red/NIR)
- Usage: once a day, 10-30cm distance

Freckle

Freckle

- Spectrum: Yellow 590nm
- Timing: 5/10/15/20 minutes
- Usage: 1-2 times a day, 5-30cm distance

Bright

- Spectrum: Red 630nm&660nm Yellow 590nm
- Timing: 5/10/15/20 minutes
- Usage: 1-2 times a day, 5-30cm distance







Shrink pores

- Spectrum: Red 630nm&660nm NIR 850nm&940nm
- Timing: 5/10/15/20 minutes
- Usage: 1-2 times a day, 5-30cm distance

Oil control & compact

- Spectrum: Blue 450nm/Yellow 590nm/Red 630nm&660nm/NIR 850nm&940nm
- Timing: 15 minutes (5 mins Blue/Yellow/Red/NIR + 10 mins Yellow/Red/NIR)
- Usage: once a day, 10-30cm distance

Settings (with glasses on):

- Long press  3 seconds to turn on the device, it will default to last used mode (first use, default to first mode).
- Press  to choose the mode, press  to set the timing, then press  to start lighting therapy, press  again to pause.
- The device will be automatically turned off when the countdown is done, and it can be turned off by long press  3 seconds manually at anytime.

Notes:

1) Make sure to consult with a physician before using the device if you have any of these conditions: acute diseases, febrile illnesses, blood pressure disorders, infectious diseases, history of skin cancer, skin disorders caused by diabetes or other diseases, currently pregnant or breast-feeding, sensitive skin due to allergic reactions, etc.

2) For better effect, please remove makeup and clean face before using the device.

3) For the initial use of red light therapy device, please set a shorter time and keep a further distance. After 3-5 times, it can be used normally.

4) Please do not look directly at the device for a long time as it may cause discomfort or irritation to the eyes.

5) The device is fragile, please use it with caution.

6) The device is not waterproof, do not use cleaning sprays. Wipe with a dry or lightly damp cloth to clean it.

Warranty clauses:

1) One year warranty from the date of purchase.

2) Free replacements are warranted for a period of three months.

3) Warranty does not cover the follow damage(s): product been modified, improper installation or use, improper clean, electrical surges, or nature caused damage.

4) Warranty does not cover the loss of third party with directly or indirect damage(s).